



PCCG Scan & Print-to-DICOM

Document integration into any current DICOM 3.0 compliant PACS System.

Print-to-DICOM Feature List

Print-to-DICOM is a powerful yet easy to use PACS utility application that converts any MS-Windows output into a DICOM 3.0 format for transmission into a PACS archive. Documents are pushed to PACS by simply selecting the configured Print-to-DICOM PACS printer.

- **MS-Windows Application Support** Converts any Windows application output into a DICOM 3.0 object for communication into any PACS product.
- **Easy to use** Select Windows PACS printer, Select and accept desired Exam and Series.
- **Simple user interface** Anyone can Print-to-DICOM in minutes without any DICOM or PACS knowledge.
- **DICOM 3.0 Compliant** Works with any modern DICOM 3.0 compliant PACS system.
- **Archive and organize presentations** Use the DICOM 3.0 architecture to organize and archive presentations.

Scan-to-DICOM Feature List

Scans paper documents and stores the pages as a series of images in your PACS system. Common documents to scan: insurance cards, lab reports, handwritten notes, etc. Using Scan-to-DICOM all relevant patient information is conveniently available at every PACS workstation for instant retrieval.

- **Scan paper documents** Any paper document can be scanned and converted into a DICOM 3.0 format for transmission and archival into a PACS architecture.
- **Acquire images from web-cameras** Pictures from any connected web-cam can be acquired and PACS archived.
- **Simple user interface** Anyone can scan documents without any DICOM or PACS knowledge.
- **Simple setup** Works with any PC networked to the PACS system.
- **Document manipulation** Allows rotation and cropping of document and images before transmission.

PRINT-to-DICOM and SCAN-to-DICOM are both DICOM Worklist Supported.

Using DICOM Worklist users eliminate human error by selecting previously established PACS systems.



PCCG, INC

Phone: (866) 279-6394
Fax: (305) 418-7433

www.pccgroup.com
sales@pccgroup.com



PCCG Scan & Print-to-DICOM

ScanToDICOM™

ScanToDICOM™ is a high-performance document scanning solution for PACS systems. ScanToDICOM™ scans any paper document and stores the pages as a series of images within the existing PACS. These paper documents - which currently impede workflow and continue the legacy of paper-based information - include imaging orders, insurance cards, lab reports, handwritten notes, etc. The solution is to use ScanToDICOM™ to have all patient-related documents and images available for viewing at every PACS workstation.

PrintToDICOM™

PrintToDICOM™ offers the ability to store computer-generated documents with the previously archived images in the PACS systems, which can then be reviewed from any PACS workstation. PrintToDICOM™ increases the PACS connectivity of any Microsoft Windows application - Microsoft Word, Excel, or Powerpoint, etc.. From Microsoft Word, for example, the selection of the Print option sends the document to the PACS instead of printing it on paper - thus avoiding continued need to file and save paper documents.

Scan-to-DICOM

Twain scan documents directly to PACS or Print digital Documents to the virtual DICOM Printer



Print-to-DICOM

Document integration into any current DICOM 3.0 compliant PACS System.



WORKLIST SUPPORT

Select or create exam and series.

A 30 day demo download may be reviewed prior to purchase.



PCCG, INC

Phone: (866) 279-6394

Fax: (305) 418-7433

www.pccgroup.com

sales@pccgroup.com